



UROLOGY

This catalog remains exclusive property of the company H. + H. Maslanka Chirurgische Instrumente GmbH.
Reprinting, even partially, is not allowed. Technical alterations and printing errors reserved.



Autoclavable at 134°C / 273°F - 5 minutes



Sterilised with gamma rays



Sterilised with ethylene oxid



Singe use only - The instruments must not be resterilised respectively reprocessed.



Radio-opaque - The instruments are visible with gamma rays.



With port for contrast medium.



With cleaning port



With high frequency connector - Please consider the generator and instrument settings.

Urology

	page
Biopsy forceps	4-11
Foreign Body Retriever	12-33
Scissors	34-41
Stone basket	42-51
Injection needle	52-56
High frequency snare	57-61
Cleaning brushes	62-66
Cytology brush	67-69
Coagulation probe	70-73
Accessory for high-frequency	74-76

Urology

Biopsy forceps



	page
Biopsy forceps - reusable - Product Group U020	5-11

Urology

Biopsy forceps

Biopsy forceps - reusable - Product Group U020

	page
Biopsy forceps with oval cups - flexible	6
Biopsy forceps with oval cups - rigid	7
Biopsy forceps with toothed cups flexible - interlocked	8
Biopsy forceps with toothed cups rigid - interlocked	9
Biopsy forceps with toothed cups flexible - successive toothed	10
Biopsy forceps with toothed cups rigid - successive toothed	11



Biopsy forceps with oval cups - flexible

- Oval cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
0,8 mm - 2,4 Fr.	on request	S0208-040	S0208-060
1,0 mm - 3 Fr.	S0210-030	S0210-040	S0210-060
1,2 mm - 3,5 Fr.	on request	S0212-040	S0212-060
1,6 mm - 5 Fr.	S0316-030	S0316-040	S0316-060
1,8 mm - 5,5 Fr.	on request	S0318-040	S0318-060
2,1 mm - 6 Fr.	on request	S0221-040	S0221-060
2,3 mm - 7 Fr.	on request	S0222-040	S0222-060

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Biopsy forceps with oval cups - rigid

- Oval cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
0,8 mm - 2,4 Fr.	on request	RS0208-040	RS0208-060
1,0 mm - 3 Fr.	RS0210-030	RS0210-040	RS0210-060
1,2 mm - 3,5 Fr.	on request	RS0212-040	RS0212-060
1,6 mm - 5 Fr.	RS0316-030	RS0316-040	RS0316-060
1,8 mm - 5,5 Fr.	on request	RS0318-040	RS0318-060
2,1 mm - 6 Fr.	on request	RS0221-040	RS0221-060
2,3 mm - 7 Fr.	on request	RS0222-040	RS0222-060

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Biopsy forceps with toothed cups flexible - interlocked

- Toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	S0410-030	S0410-040	S0410-060
1,6 mm - 5 Fr.	S0416-030	S0416-040	S0416-060
1,8 mm - 5,5 Fr.	S0418-030	S0418-040	S0418-060
2,1 mm - 6,3 Fr.	S0421-030	S0421-040	S0421-060
2,3 mm - 7 Fr.	S0422-030	S0422-040	S0422-060

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Biopsy forceps with toothed cups rigid - interlocked

- Toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	RS0410-030	RS0410-040	RS0410-060
1,6 mm - 5 Fr.	RS0416-030	RS0416-040	RS0416-060
1,8 mm - 5,5 Fr.	RS0418-030	RS0418-040	RS0418-060
2,1 mm - 6,3 Fr.	RS0421-030	RS0421-040	RS0421-060
2,3 mm - 7 Fr.	RS0422-030	RS0422-040	RS0422-060

Required accessories: Handle and cable
Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Biopsy forceps with toothed cups flexible - successive toothed

- Toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,6 mm - 5 Fr.	S04416-030	S04416-040	S04416-060
2,1 mm - 6,3 Fr.	S04421-030	S04421-040	S04421-060
2,3 mm - 7 Fr.	S04422-030	S04422-040	S04422-060

Required accessories: Handle and cable
Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Biopsy forceps with toothed cups rigid - successive toothed

- Toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,6 mm - 5 Fr.	RS04416-030	RS04416-040	RS04416-060
2,1 mm - 6,3 Fr.	RS04421-030	RS04421-040	RS04421-060
2,3 mm - 7 Fr.	RS04422-030	RS04422-040	RS04422-060

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"

Urology

Foreign Body Retriever




	page
Foreign Body Retriever - reusable - Product Group U022	13-33

Urology

Foreign Body Retriever

Foreign Body Retriever - reusable - Product Group U022

	page
Foreign Body Retriever serrated jaws - 4,8 mm flexible - successive toothed	14
Foreign Body Retriever serrated jaws - 4,8 mm rigid - successive toothed	15
Foreign Body Retriever serrated jaws 6,6 mm long - flexible - interlocked	16
Foreign Body Retriever serrated jaws 6,6 mm - rigid - interlocked toothed	17
Foreign Body Retriever - flexible - interlocked toothed	18
Foreign Body Retriever - rigid - interlocked toothed	19
Foreign Body Retriever - one sided flexible - interlocked toothed	20
Foreign Body Retriever - one sided rigid - interlocked toothed	21
Foreign Body Retriever - flexible	22
Foreign Body Retriever - rigid	23
Foreign Body Retriever for stone-street - flexible	24
Foreign Body Retriever for stone-street - rigid	25
DJ - Forceps - flexible	26
DJ - Forceps - rigid	27
Rat tooth forceps - flexible	28
Rat tooth forceps - rigid	29
Foreign Body Retriever - flexible	30
Foreign Body Retriever - rigid	31
Forceps for stones and fragments - rigid	32
Forceps for stones and fragments - type curved	33




**Foreign Body Retriever serrated jaws -
4,8 mm flexible - successive toothed**

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,2 mm - 3,5 Fr.	on request	S04612-040	S04612-060
1,6 mm - 5 Fr.	on request	S04616-040	S04616-060



Foreign Body Retriever serrated jaws - 4,8 mm rigid - successive toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,2 mm - 3,5 Fr.	on request	RS04612-040	RS04612-060
1,6 mm - 5 Fr.	on request	RS04616-040	RS04616-060



Foreign Body Retriever serrated jaws 6,6 mm long - flexible - interlocked

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,6 mm - 5 Fr.	on request	R04316-040	R04316-060
2,1 mm - 6 Fr.	on request	R04321-040	R04321-060




Foreign Body Retriever serrated jaws 6,6 mm - rigid - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,6 mm - 5 Fr.	on request	RS04316-040	RS04316-060
2,1 mm - 6 Fr.	on request	RS04321-040	RS04321-060




Foreign Body Retriever - flexible - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,2 mm - 3,5 Fr.	on request	S05212-040	S05212-060
2,3 mm - 7 Fr.	on request	S05222-040	S05222-060



Foreign Body Retriever - rigid - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,2 mm - 3,5 Fr.	on request	RS05212-040	RS05212-060
2,3 mm - 7 Fr.	on request	RS05222-040	RS05222-060



Foreign Body Retriever - one sided flexible - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
	40 cm
2,3 mm - 7 Fr.	S05322-040



Foreign Body Retriever - one sided rigid - interlocked toothed

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
	40 cm
2,3 mm	RS05322-040



Foreign Body Retriever - flexible

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	S1810-030	S1810-040	S1810-060
1,2 mm - 3,5 Fr.	on request	S1812-040	S1812-060
1,6 mm - 5 Fr.	S1816-030	S1816-040	S1816-060
1,8 mm - 5,5 Fr.	on request	S1818-040	S1818-060
2,4 mm - 7 Fr.	S1822-030	S1822-040	S1822-060

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Foreign Body Retriever - rigid

- Cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	RS1810-030	RS1810-040	RS1810-060
1,2 mm - 3,5 Fr.	on request	RS1812-040	RS1812-060
1,6 mm - 5 Fr.	RS1816-030	RS1816-040	RS1816-060
1,8 mm - 5,5 Fr.	on request	RS1818-040	RS1818-060
2,4 mm - 7 Fr.	RS1822-030	RS1822-040	RS1822-060

Required accessories: Handle and cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Foreign Body Retriever for stone-street - flexible

- Long cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
1,6 mm - 5 Fr.	60 cm S04116-060




Foreign Body Retriever for stone-street - rigid

- Long cross toothed cups
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
1,6 mm - 5 Fr.	60 cm RS04116-060



DJ - Forceps - flexible

- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	40 cm	60 cm
2,3 mm - 7 Fr.	S08322-040	S08322-060
3,0 mm - 9 Fr.	S08330-040	S08330-060



DJ - Forceps - rigid

- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	40 cm	60 cm
2,3 mm - 7 Fr.	RS08322-040	RS08322-060
3,0 mm - 9 Fr.	RS08330-040	RS08330-060



Rat tooth forceps - flexible

- 2:1 teeth
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	40 cm	60 cm
1,6 mm - 5 Fr.	S08416-040	S08416-060
2,3 mm - 7 Fr.	S08422-040	S08422-060



Rat tooth forceps - rigid

- 2:1 teeth
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	40 cm	60 cm
1,6 mm - 5 Fr.	RS08416-040	RS08416-060
2,3 mm - 7 Fr.	RS08422-040	RS08422-060

Foreign Body Retriever - flexible



- Forked jaws with 1:2 tooth
- Jaws not totally retractable
- Handle with cleaning port



Instrument diameter	Working length
	60 cm
1,6 mm - 5 Fr.	S1016-060
1,8 mm - 5,5 Fr.	S1018-060


Foreign Body Retriever - rigid



- Forked jaws with 1:2 tooth
- Jaws not totally retractable
- Handle with cleaning port



Instrument diameter	Working length
1,6 mm - 5 Fr.	60 cm RS1016-060
1,8 mm - 5,5 Fr.	RS1018-060



Forceps for stones and fragments - rigid

- Fenestrated
- Ground spiral
- handle with cleaning port



Instrument diameter	Working length
3,0 mm - 9 Fr.	40 cm RS05130-040



**Forceps for stones and fragments -
type curved**

- Type curved - rigid
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
	30 cm
3,0 mm - 9 Fr.	RS18130-030

Urology

Scissors



	page
Scissors - reusable - Product Group U024	35-41

Urology

Scissors

Scissors - reusable - Product Group U024

	page
Scissors one-sided - blunt tips flexible	36
Scissors one-sided - blunt tips rigid	37
Scissors double-sided - blunt tips rigid	38
Scissors one-sided - pointed tips flexible	39
Scissors one-sided - pointed rigid	40
Hook Scissors, single action	41



Scissors one-sided - blunt tips flexible

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	S09103-030	S09103-040	S09103-060
1,6 mm - 5 Fr.	S09163-030	S09163-040	S09163-060
2,1 mm - 6 Fr.	S09213-030	S09213-040	S09213-060
2,3 mm - 7 Fr.	S09223-030	S09223-040	S09223-060



Scissors one-sided - blunt tips rigid

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length		
	30 cm	40 cm	60 cm
1,0 mm - 3 Fr.	RS09103-030	RS09103-040	RS09103-060
1,6 mm - 5 Fr.	RS09163-030	RS09163-040	RS09163-060
2,1 mm - 6 Fr.	RS09213-030	RS09213-040	RS09213-060
2,3 mm - 7 Fr.	RS09223-030	RS09223-040	RS09223-060



**Scissors double-sided - blunt tips
rigid**

- Double-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	30 cm	40 cm
1,6 mm - 5 Fr.	RS0916-030	RS0916-040



Scissors one-sided - pointed tips flexible

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port
- Pointed tips



Instrument diameter	Working length	
	30 cm	40 cm
1,6 mm - 5 Fr.	S09162-030	S09162-040
2,3 mm - 7 Fr.	S09222-030	S09222-040



Scissors one-sided - pointed rigid

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length	
	30 cm	40 cm
1,6 mm - 5 Fr.	RS09162-030	RS09162-040
2,3 mm - 7 Fr.	RS09222-030	RS09222-040



Hook Scissors, single action

- One-sided movable scissors
- Ground spiral
- Handle with cleaning port



Instrument diameter	Working length
2,3 mm - 7 Fr.	70 cm
	RS09224

Urology

Stone basket



page

Stone basket - single use - Product Group U031	43-45
Stone basket - reusable - Product Group U030	46-51

Urology

Stone basket

Stone basket - single use - Product Group U031

	page
Stone basket 4-fold turned - single use Ø 1,65 mm - 5 Fr.	44
Stone basket 4-fold turned - single use Ø 1,8 mm - 5,5 Fr.	45



Stone basket 4-fold turned - single use

Ø 1,65 mm - 5 Fr.

For working channel Ø 1,8 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Fixed handle



Basket type	Basket size	Working length
		60 cm
Monofile wire	40 mm	E150164-060
	50 mm	E150165-060

Also available in nitinol! Please enquire.



Stone basket 4-fold turned - single use

Ø 1,8 mm - 5,5 Fr.

For working channel Ø 2,0 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Fixed handle



Basket type	Basket size	Working length
		60 cm
Monofile wire	40 mm	E150184-060
	50 mm	E150185-060

Also available in nitinol! Please enquire.

Urology

Stone basket

Stone basket - reusable - Product Group U030

	page
Stone basket 4-fold turned Ø 1,0 mm - 3 Fr.	47
Stone basket 4-fold turned Ø 1,65 mm - 5 Fr.	48
Stone basket 4-fold turned Ø 1,8 mm - 5,5 Fr.	49
Stone basket 6-fold turned Ø 1,65 mm - 5 Fr.	50
Stone basket 6-fold turned Ø 1,8 mm - 5,5 Fr.	51



Stone basket 4-fold turned

Ø 1,0 mm - 3 Fr.

For working channel Ø 1,2 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket Type	Basket size	Working length	
		40 cm	60 cm
Monofile wire	40 mm	150104-040	150104-060



Stone basket 4-fold turned

Ø 1,65 mm - 5 Fr.

For working channel Ø 1,8 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket type	Basket size	Working length	
		40 cm	60 cm
Monofile wire	40 mm	150164-040	150164-060
	50 mm	150165-040	150165-060
Multifile wire	40 mm	160164-040	160164-060
	50 mm	160165-040	160165-060



Stone basket 4-fold turned Ø 1,8 mm - 5,5 Fr.

For working channel Ø 2,0 mm and higher

- Twisted construction
- 4 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket type	Basket size	Working length	
		40 cm	60 cm
Monofile wire	40 mm	150184-040	150184-060
	50 mm	150185-040	150185-060
Multifile wire	40 mm	160184-040	160184-060
	50 mm	160185-040	160185-060



Stone basket 6-fold turned Ø 1,65 mm - 5 Fr.

For working channel Ø 1,8 mm and higher

- Twisted construction
- 6 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket type	Basket size	Working length	
		40 cm	60 cm
Monofile wire	40 mm	152164-040	152164-060
	50 mm	152165-040	152165-060
Multifile wire	40 mm	162164-040	162164-060
	50 mm	162165-040	162165-060



Stone basket 6-fold turned

Ø 1,8 mm - 5,5 Fr.

For working channel Ø 2,0 mm and higher


- Twisted construction
- 6 wires
- With Luer Lok port for contrast medium and cleaning
- Can be taken apart for cleaning
- With handle



Basket type	Basket size	Working length	
		40 cm	60 cm
Monofile wire	40 mm	152184-040	152184-060
	50 mm	152185-040	152185-060
Multifile wire	40 mm	162184-040	162184-060
	50 mm	162185-040	162185-060

Urology

Injection needle



	page
Injection needle - single use - Product Group U001	53-54
Injection needle - reusable - Product Group U006	55-56

Urology

Injection needle

Injection needle - single use - Product Group U001

	page
Injection needle with needle - 5,5 Fr. single use	54



Injection needle with needle - 5,5 Fr. single use

For working channels Ø 2,0 mm and higher

- Needle diameter 0,7 mm
- Needle length 4 mm or 6 mm
- Sterile packed



Needle length	Working Length	Type
		Plastic tube
4 mm	60 cm	ET2518.74-060
6 mm	60 cm	ET2518.76-060

Urology

Injection needle

Injection needle - reusable - Product Group U006

	page
Injection needle with needle - 5,5 Fr.	56



Injection needle with needle - 5,5 Fr.


- Needle diameter 0,7 mm
- 60 cm working length



Needle length	Model	Type
		Plastic tube
10 mm	Needle retractable	2518-50/10
10 mm	Fixed	2518-50/10F

Urology

High frequency snare



page


High frequency snare - single use - Product Group U00558-61

Urology

High frequency snare

High frequency snare - single use - Product Group U005

	page
Symmetrical High Frequency Snare Ø 1,65 mm - 5 Fr. - single use	59
Symmetrical High Frequency Snare Ø 1,8 mm - 5,5Fr. - single use	60
Symmetrical High Frequency Snare Ø 2,35 mm - 7 Fr. - single use	61



Symmetrical High Frequency Snare Ø 1,65 mm - 5 Fr. - single use

For working channels Ø 1,8 mm and higher


- Symmetrical snare-opening
- Available in two loop sizes
- With Handle



HF connection type	Working length	Loop size	
		20 mm	25 mm
4 mm socket	60 cm	E341615-060	E341616-060
Olympus	60 cm	E341615A-060	E341616A-060

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Symmetrical High Frequency Snare Ø 1,8 mm - 5,5Fr. - single use

For working channels Ø 2,0 mm and higher

- Symmetrical snare-opening
- Available in two loop sizes
- With Handle



HF connection type	Working length	Loop size	
		20 mm	25 mm
4 mm socket	60 cm	E341815-060	E341816-060
Olympus	60 cm	E341815A-060	E341816A-060

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Symmetrical High Frequency Snare Ø 2,35 mm - 7 Fr. - single use

For working channels Ø 2,45 mm and higher

- Symmetrical snare-opening
- Available in two loop sizes
- With Handle



HF connection type	Working length	Loop size	
		20 mm	25 mm
4 mm socket	60 cm	E342215-060	E342216-060
Olympus	60 cm	E342215A-060	E342216A-060

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"

Urology

Cleaning brushes



	page
Cleaning brush single use - Product Group U004	63-64
Cleaning brush - reusable - Product Group U004	65-66


Urology

Cleaning brushes

Cleaning brush single use - Product Group U004

	page
Cleaning brush - single use	64

Cleaning brush - single use



- With double-sided brush
- For different working channels



	Working length	For working channel	
		Ø 1,65 mm and higher	Ø 2,0 mm and higher
Double-sided	100 cm	E4116-100	E4118-100

Urology

Cleaning brushes

Cleaning brush - reusable - Product Group U004

	page
Cleaning brush	66

Cleaning brush



- With single oder double-sided brush
- Flexible distal tip
- For different working channels



	Working length	For working channel			
		Ø 1,2 mm and higher	Ø 1,65 mm and higher	Ø 2,0 mm and higher	Ø 2,6 mm and higher
Single-sided	40 cm	41101-040	41161-040	41181-040	41221-040
	60 cm	41101-060	41161-060	41181-060	41221-060
Double-sided	40 cm	BS41101-040	BS41161-040	BS41181-040	BS41221-040
	60 cm	BS41101-060	BS41161-060	BS41181-060	BS41221-060

Urology

Cytology brush



	page
Cytology brush for single use - Product Group U003	68-69

Urology

Cytology brush

Cytology brush for single use - Product Group U003

	page
Cytology brush for single-use	69



Cytology brush for single-use

- With plastic tube
- Brush totally retractable
- With ball for tube sealing
- Single sterile packed
- Packing unit 10 pieces



Working length	For working channel	
	Ø 1,8 mm and higher	Ø 2,0 mm and higher
120 cm	E4216-120	E4218-120

Urology

Coagulation probe



	page
Coagulation probe -reusable - Product Group U007	71-73

Urology

Coagulation probe

Coagulation probe -reusable - Product Group U007

	page
Coagulation probe Ø 1,65 mm - 5 Fr.	72
Coagulation probe Ø 1,8 mm - 5,5 Fr.	73



Coagulation probe Ø 1,65 mm - 5 Fr.

For working channels Ø 1,8 mm and higher

- With or without Axial flushing port
- Flexible



Spiral type	Working length	Art. No.
With flushing port	40 cm	8016-040
	60 cm	8016-060
Without flushing port	40 cm	80161-040
	60 cm	80161-060

Required accessories: Cable

Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"



Coagulation probe Ø 1,8 mm - 5,5 Fr.

For working channels Ø 2,0 mm and higher

- With or without Axial flushing port
- Flexible



Spiral type	Working length	Art.No.
With flushing port	40 cm	8018-040
	60 cm	8018-060
Without flushing port	40 cm	80181-040
	60 cm	80181-060

Required accessories: Cable
Please find the appropriate article numbers on our accessory pages "accessory for high-frequency"

Urology

Accessory for high-frequency



	page
High frequency cable - Product Group 0033	75-76

Urology

Accessory for high-frequency

High frequency cable - Product Group 0033

	page
High frequency cable	76

High frequency cable



- 4 mm connector (Martin, Erbe, Storz)
- 5 mm connector (Erbe)
- 6 mm connector (Olympus)
- 8 mm connector (Vallylab)



Generator side	Handle side	
	4 mm socket	Olympus
Martin Erbe (ACC) Storz	6400	6400A
Erbe (ICC)	6404	6404A
Olympus	6402	6402A
Vallylab	6407	6407A